

ARM SELECTIONS Models 2615, 2614, 2601, 2605



Rail Mount Accessories

slide along BOYD's "On Track Rail System". The Rail Mech allows front to back, up and down, and rotational positioning. Loosen the knob and position, tighten knob and you are locked into place. All Rail Mounted Accessories are universal and mount right or left. 2615, 2614 Models are equipped standard with Seat Rails. Back Rails are OPTIONAL. 2601, 2605 Models are equipped standard with Seat and Back Rails.



PA-110 Patient Arm
Universal Mount Right or
Left on BOYD's "On Track
Rail System"
P/N 68-6223 Putty Trim
P/N 68-6223B Black Trim




IV-119 Multi Position IV Arm
Board Universal Mount
Right or Left on BOYD's "On
Track Rail System"
P/N 68-6266 Putty Trim
P/N 68-6266B Black Trim



IV-123 Multi Position IV Arm
Board Universal Mount Right
or Left on BOYD's "On Track
Rail System"
P/N 68-6271 Putty Trim
P/N 68-6271B Black Trim

ARM SELECTIONS Models 2615, 2614



Compensating Arms

automatically adjust with chair back. All arms are interchangeable from side to side, removable, and width adjustable. IV Arms rotate in multi positions with the touch of a lever. Arms drop to the side and out of the way for patient access/egress. A firm locking action holds the arm stable during surgical procedures.



PA-220 Patient Arm;
COMPENSATING
Universal Mount Right or Left
into BOYD Back Hub
P/N 68-6224 Putty Trim
P/N 68-6224B Black Trim



IV-219 Multi Position IV Arm
Board; COMPENSATING
Universal Mount Right or Left
into BOYD Back Hub
P/N 68-6272 Putty Trim
P/N 68-6272B Black Trim